

20X LumiGLO[®] Reagent and 20X Peroxide

- ✓ 25 ml 20X LumiGLO (luminol chemiluminescent substrate) and 25 ml 20X Peroxide

Orders ■ 877-616-CELL (2355)
orders@cellsignal.com


Support ■ 877-678-TECH (8324)
info@cellsignal.com

Web ■ www.cellsignal.com

rev. 11/09/09

This product is for *in vitro* research use only and is not intended for use in humans or animals. This product is not intended for use as a therapeutic or in diagnostic procedures.

Background: Chemiluminescent detection systems have emerged as the best all-around method for detection of Western blots. They eliminate the hazards associated with radioactive materials and toxic chromogenic substrates. The speed and sensitivity of these methods are unequalled by traditional alternatives. Because results are generated on film, it is possible to record and store data permanently, and blots detected with chemiluminescent methods are easily stripped for subsequent reprobing with additional antibodies. Horseradish peroxidase (HRP) conjugated secondary antibodies are utilized in conjunction with specific chemiluminescent substrates to generate the light signal. Horseradish peroxidase-antibody conjugates have a very high turnover rate, giving good sensitivity with short reaction times.



After the primary antibody is bound to the target protein, a complex with HRP-linked secondary antibody is formed. The LumiGLO[®] is added and emits light during enzyme catalyzed decomposition.

Description: LumiGLO[®] chemiluminescent substrate is a luminol-based system designed for use with our Phototope[®]-HRP detection assays utilizing peroxidase-labeled antibodies immobilized on membranes. In the presence of hydrogen peroxide, horseradish peroxidase (HRP) converts luminol to an excited intermediate dianion. This dianion emits light on return to its ground state. Light emission is maximal immediately after exposure of the substrate to HRP and continues for 0.5-1 hour. Light can be captured on X-ray film, typically by exposure for a few seconds. Maximum sensitivity can be obtained by longer exposure.

*Avoid repeated exposure to skin (see enclosed Material Safety Data Sheet or refer to our website for further information).

Direction for Use:

- Wash membrane-bound HRP (antibody conjugate) three times, for 5 minutes in TBS/T.
- Prepare substrate by diluting 20X LumiGLO[®] and 20X Peroxide to 1X in water (e.g. For 10 ml, add 0.5 ml LumiGLO[®] and 0.5 ml peroxide to 9.0 ml water).
- Incubate substrate with membrane for 1 minute, remove excess solution (membrane remains wet), wrap in plastic and expose to X-ray film.

Solutions and Reagents:

Note: Prepare solutions with Milli-Q[®] or equivalently purified water.

Wash Buffer (TBS/T): 20 mM Tris-HCl (pH 7.6), 137 mM NaCl and 0.1% Tween-20

Advantages of our Phototope[®] Western Detection Systems
-Sensitivity: Detection of sub- picogram amounts of protein

Reaction Conditions:

- Wash membrane-bound HRP (antibody conjugate) three times for 5 minutes in TBS/T.
- Prepare substrate by diluting 20X LumiGLO and 20X Peroxide to 1X in water (e.g., for 10 ml, add 0.5 ml LumiGLO and 0.5 ml peroxide to 9.0 ml water).
- Incubate substrate with membrane for 1 minute, remove excess solution (membrane remains wet), wrap in plastic and expose to X-ray film.

Storage: Store at 4°C.

Companion Products:

Phototope[®]-HRP Western Blot Detection System, Anti-rabbit IgG, HRP-linked Antibody #7071

Phototope[®]-HRP Western Blot Detection System, Anti-mouse IgG, HRP-linked Antibody #7072

Anti-rabbit IgG, HRP-linked Antibody #7074

Anti-biotin, HRP-linked Antibody #7075

Anti-mouse IgG, HRP-linked Antibody #7076

Anti-rat IgG, HRP-linked Antibody #7077

Advantages of our Phototope[®] Western Detection Systems

• **Sensitivity:** Detection of subpicogram amounts of protein is routine with good primary antisera.

• **Speed:** Less than 1 hour is required for the entire detection procedure. Exposure times are seconds to minutes.

• **Multiple Exposures:** Light is emitted at a constant rate for several minutes, so you can perform multiple exposures to optimize signal intensity. Re-exposure at a future date is achieved by simply adding more reagent.

• **Stability:** A permanent hard-copy record is generated that will not fade or disintegrate over time.

• **Quantitative:** X-ray films can be scanned to quantitate band intensities.

*Avoid repeated exposure to skin (see enclosed Material Safety Data Sheet or refer to our website for further information).

Phototope[®] is a trademark of Cell Signaling Technology, Inc. LumiGLO[™] is a trademark of Kirkegaard & Perry Laboratories (KPL). Milli-Q[™] is a trademark of Millipore.

Material Safety Data Sheet (MSDS) for 20X LumiGLO® and 20X Peroxide

I. Identification:

Product name: 20X LumiGLO® and 20X Peroxide
Product Catalog: 7003
CAS number: None
Manufacturer Supplier: Cell Signaling Technology
 3 Trask Lane
 Danvers, MA 01923 USA
 1-978-867-2300 TEL
 1-978-867-2400 FAX
 1-978-578-6737 Emergency TEL

II. Composition/Information on Ingredients:

| Hazardous Reagent: | Percent | CAS# |
|--------------------|---------|---------|
| Dimethylsulfoxide | ≤20% | 67-68-5 |

This product is For Research Use Only. According to 29 CFR 1910.1200(d), mixtures with hazardous ingredients at less than <1% and carcinogens at less than <0.1% are considered non-hazardous.

III. Hazard Identification:

CAUTION: This product is not for use in humans. It is intended for research purposes only. To the best of our knowledge, the chemical, physical, and toxicological properties of this material have not been established.

Emergency Overview: Irritant. Irritating to eyes, respiratory system, skin.

Potential Health Effects:

Inhalation: May be harmful if inhaled. May cause respiratory tract irritation.

Eye Contact: May cause eye irritation.

Skin Contact: May be harmful if absorbed through skin. Prolonged or repeated contact may cause skin irritation.

Ingestion: May be harmful if swallowed.

IV. First Aid Measures:

Inhalation: Remove to fresh air. If breathing is difficult, get medical attention.

Ingestion: If person is conscious, wash out mouth with water. Get medical attention.

Skin exposure: Wash skin with soap and water. If irritation develops or persists, get medical attention.

Eye exposure: Immediately flush eyes water for at least 15 minutes. Get medical attention.

V. Fire Fighting Measures:

Flash Point: N/A

Autoignition Temperature: N/A

Explosion: N/A

Fire extinguishing media: water spray, dry chemical, alcohol foam, or carbon dioxide.

Firefighting: wear protective clothing and self-contained breathing apparatus to prevent contact with skin and eyes. May emit toxic fumes under fire conditions.

VI. Accidental Release Measures:

Wear appropriate personal protective equipment as indicated in Section VIII. Absorb liquid with an absorbent material. Transfer contaminated absorbent to a closed chemical waste container for disposal. Wash spill site after material has been picked up for disposal.

VII. Handling And Storage:

Store at 4°C in tightly closed container.

Avoid inhalation of vapor or mist. Avoid contact with eyes, skin, and clothing. Wash thoroughly after handling. Avoid prolonged or repeated exposure.

VIII. Exposure Controls/Personal:

Ventilation System: a system of local and/or general exhaust is recommended.

Skin Protection: wear compatible chemical resistant gloves and protective clothing.

Eye protection: wear protective safety glasses or chemical safety goggles. Maintain eye wash fountain and quick-drench facilities in work area.

IX. Physical And Chemical Properties:

| | |
|--|-----------------------------------|
| Appearance: | clear faint yellow colored liquid |
| Odor: | data not available |
| pH: | data not available |
| Boiling Point: | >100°C/212°F (water) |
| Melting or Freezing Point: | <0°C/32°F (water) |
| Flash Point: | data not available |
| Volatile Organic Compounds (VOC): | data not available |
| Autoignition temp. | data not available |
| Solubility (water) | miscible in water |

X. Stability and Reactivity:

Stability: Stable under normal conditions.

Conditions to avoid: strong oxidizing agents, strong acids, strong bases.

Hazardous Decomposition: carbon monoxide, carbon dioxide.

Hazardous polymerization: will not occur.

XI. Toxicological Information:

Acute toxicity: data not available. Chronic exposure: data not available

Potential Health Effects:

Inhalation: May be harmful if inhaled. Causes respiratory tract irritation.

Skin: May be harmful if absorbed through skin. Causes skin irritation.

Eyes: Causes eye irritation.

Ingestion: Harmful if swallowed.

Toxicity Data on Hazardous ingredient Dimethyl Sulfoxide, CAS#67-68-5

RTECS: PV6210000

LD50 Oral - rat - 14,500 mg/kg

LC50 Inhalation - rat - 4 h - 40250 ppm

LD50 Dermal - rabbit - > 5,000 mg/kg

XII. Ecological Information:

No data available.

XIII. Disposal Considerations: Dispose of in accordance with federal, state and local environmental regulations. Contact a licensed professional waste disposal service to dispose of this material.

XIV. Transport Information:

D.O.T. Proper Shipping Name: None. This substance is considered non-hazardous for transport.

IATA Proper Shipping Name: None. This substance is considered non-hazardous for air transport.

XV. Regulatory Information:

EU: Not classified

OSHA: Ingredient Dimethyl Sulfoxide, CAS#67-68-5: Combustible Liquid, Target Organ Effect

Canadian DSL: Listed: Ingredient Dimethyl Sulfoxide, CAS#67-68-5

SARA 302, 313 Ingredients Not Listed.

SARA 311/312: Ingredient Dimethyl Sulfoxide, CAS#67-68-5: Fire Hazard, Chronic Health Hazard.

Massachusetts Right To Know: Ingredients Not Listed.

Pennsylvania Right To Know: Ingredient Dimethyl Sulfoxide, CAS#67-68-5

New Jersey Right To Know: Ingredient Dimethyl Sulfoxide, CAS#67-68-5

California Prop. 65: Ingredients Not Listed.

XVI. Other Information:

This product is for research use only and is not intended for use in humans. To the best of our knowledge, this document is accurate. It is intended to serve as a guide for safe use of this product in a laboratory setting by experienced personnel. The burden of safe use of this material rests entirely with the user. The above information is believed to be accurate but is not necessarily all-inclusive and shall be used only as a guide. Cell Signaling Technology, Inc., shall not be held liable for any damage resulting from the handling of or from contact with the above product.