

PathScan® Phospho-Akt (Thr308) Sandwich ELISA kit

✓ 1 Kit
(96 assays)

Orders ■ 877-616-CELL (2355)
orders@cellsignal.com

Support ■ 877-678-TECH (8324)
info@cellsignal.com

Web ■ www.cellsignal.com

rev. 05/29/09

This product is for *in vitro* research use only and is not intended for use in humans or animals.
This product is not intended for use as a therapeutic or in diagnostic procedures.

Species Cross-Reactivity: H, M

Description: CST's PathScan® Phospho-Akt (Thr308) Sandwich ELISA Kit is a solid phase sandwich enzyme-linked immunosorbent assay (ELISA) that detects phospho-Akt (Thr308) protein. An Akt Antibody has been coated onto the microwells. After incubation with cell lysates, both phospho- and nonphospho-Akt proteins are captured by the coated antibody. Following extensive washing, Phospho-Akt (Thr308) Mouse mAb is added to detect the captured phospho-Akt protein. Anti-mouse IgG, HRP-linked Antibody is then used to recognize the bound detection antibody. HRP substrate, TMB, is added to develop color. The magnitude of optical density for this developed color is proportional to the quantity of phospho-Akt (Thr308) protein.

Please visit www.cellsignal.com for a complete listing of recommended companion products.

Specificity/Sensitivity: CST's PathScan® Phospho-Akt (Thr308) Sandwich ELISA Kit #7252 detects endogenous levels of phospho-Akt (Thr308) protein. As shown in Figure 1, using Phospho-Akt (Thr308) ELISA Kit #7252, a significant induction of Phospho-Akt (Thr308) is detected in NIH/3T3 cells treated with PDGF. However, levels of total Akt (phospho and nonphospho) detected by PathScan® Total Akt Sandwich ELISA Kit #7170, remain unchanged (Figure 1). Phospho-Akt (Thr308) in Jurkat cells is also detected by this ELISA kit #7252.

Products Included	Volume	Solution Color
Akt Antibody Coated Microwells*	96 tests	
Phospho-Akt (Thr308) Detection Ab	11 ml	green
Anti-mouse IgG HRP-Linked Ab	11 ml	red
TMB Substrate	11 ml	colorless
STOP Solution	11 ml	colorless
Sealing Tape	2 sheets	
20X Wash Buffer	25 ml	colorless
Sample Diluent	25 ml	blue
10X Cell Lysis Buffer #9803**	15 ml	yellowish

* 12 8-well modules -Each module is designed to break apart for 8 tests.

**Kit should be stored at 4°C with the exception of 10X Cell Lysis Buffer, which is stored at -20°C (packaged separately).




Figure 1: Treatment of NIH/3T3 cells with PDGF stimulates phosphorylation of Akt at Thr308, detected by PathScan® Phospho-Akt (Thr308) Sandwich ELISA Kit #7252, but does not affect the level of total Akt1 detected by PathScan® Total Akt1 Sandwich ELISA Kit #7170. OD 450 nm readings are shown in the top figure, while the corresponding Western blots using Phospho-Akt (Thr308) Antibody #9275 (right panel) or Akt Antibody #9272 (left panel), are shown in the bottom figure.



Background: Akt, also referred to as PKB or Rac, plays a critical role in controlling survival and apoptosis (1-3). This protein kinase is activated by insulin and various growth and survival factors and functions in a wortmannin-sensitive pathway involving PI3 kinase (2,3). Akt is activated by phospholipid binding and activation loop phosphorylation at Thr308 by PDK1 (4) and by phosphorylation within the carboxy terminus at Ser473. The previously elusive PDK2 responsible for phosphorylation of Akt at Ser473 has been identified as mammalian target of rapamycin (mTor) in a rapamycin-insensitive complex with Rictor and Sin1, termed TORC2 (5,6). Akt promotes cell survival by inhibiting apoptosis through its ability to phosphorylate and inactivate several targets, including Bad (7), Forkhead transcription factors (8), c-Raf (9) and caspase-9. PTEN phosphatase is a major negative regulator of the PI3 kinase/Akt signaling pathway (10). LY294002 is a specific PI3 kinase inhibitor (11).

One of the essential functions of Akt is the regulation of glycogen synthesis through phosphorylation and inactivation of GSK-3 α and β (12,13). Akt may also play a role in insulin stimulation of glucose transport (12).

In addition to its role in survival and glycogen synthesis, Akt is involved in cell cycle regulation by preventing GSK-3 β mediated phosphorylation and degradation of cyclin D1 (14) and by negatively regulating the cyclin dependent kinase inhibitors p27 Kip (15) and p21 Waf1 (16). Akt also plays a critical role in cell growth by directly phosphorylating mTOR in a rapamycin-sensitive complex containing raptor, termed TORC1 (17). More importantly, Akt phosphorylates and inactivates tuberlin (TSC2), an inhibitor of mTOR within the TORC1 complex (18). Inhibition of mTOR stops the protein synthesis machinery due to inactivation of its effector, p70 S6 kinase and activation of the eukaryotic initiation factor, 4E binding protein 1 (4E-EP1), an inhibitor of translation (19,20).

Background References:

- (1) Franke, T.F. (1997) *Cell* 88, 435–437.
- (2) Burgering, B.T. and Coffey, P.J. (1995) *Nature* 376, 599–602.
- (3) Franke, T.F. et al. (1995) *Cell* 81, 727–736.
- (4) Alessi, D.R. et al. (1996) *EMBO J.* 15, 6541–6551.
- (5) Sarbassov, D.D. et al. (2005) *Science* 307, 1098–1101.
- (6) Jacinto, E. et al. (2006) *Cell* 127, 125–137.
- (7) Cardone, M.H. et al. (1998) *Science* 282, 1318–1321.
- (8) Brunet, A. et al. (1999) *Cell* 96, 857–868.
- (9) Zimmerman, S. et al. (1999) *Science* 286, 1741–1744.
- (10) Cantley, L.C. et al. (1999) *Proc. Natl. Acad. Sci. USA* 96, 4240–4245.
- (11) Vlahos, C. et al. (1994) *J. Biol. Chem.* 269, 5241–5248.
- (12) Hajdich, E. et al. (2000) *FEBS Lett.* 492, 199–203.
- (13) Cross, D.A. et al. (1995) *Nature* 373, 785–789.
- (14) Diehl, J.A. et al. (1998) *Genes Dev.* 12, 3499–3511.
- (15) Gesbert, F. et al. (2000) *J. Biol. Chem.* 275, 39223–39230.
- (16) Zhou, B.P. et al. (2001) *Nat. Cell Biol.* 3, 245–252.
- (17) Nave, B.T. et al. (1999) *Biochem. J.* 344, 427–431.
- (18) Manning, B.D. et al. (2000) *Mol. Cell* 4, 648–657.
- (19) Manning, B.D. et al. (2002) *Mol. Cell* 10, 151–162.
- (20) Inoki, K. et al. (2002) *Nat. Cell Biol.* 4, 648–657.




Figure 2: The relationship between protein concentration of lysates from untreated and PDGF-treated NIH/3T3 cells and kit assay optical density readings is shown. After starvation, NIH/3T3 cells (85% confluence) were treated with PDGF (50 ng/ml) for 10 min at 37°C, and then lysed.

Sandwich ELISA Protocol

A Reagent Preparation

1. Bring all microwell strips to room temperature before use.
2. Prepare 1X Wash Buffer by diluting 20X Wash Buffer (included in each PathScan® Sandwich ELISA Kit) in Milli-Q or equivalently purified water.
3. **1X Cell Lysis Buffer from CST #9803:** 20 mM Tris (pH 7.5), 150 mM NaCl, 1 mM ethylene diamine tetraacetate (EDTA), 1 mM ethylene glycol-bis(2-aminoethyl)-N,N,N',N'-tetraacetic acid (EGTA), 1% Triton X-100, 2.5 mM sodium pyrophosphate, 1 mM β-glycerophosphate, 1 mM Na₃VO₄, 1 μg/ml leupeptin. This buffer can be stored at 4°C for short-term use (1–2 weeks).

B Preparing Cell Lysates

1. Aspirate media. Treat cells by adding fresh media containing regulator for desired time.
2. To harvest cells under nondenaturing conditions, remove media and rinse cells once with ice-cold PBS.
3. Remove PBS and add 0.5 ml ice-cold 1X Cell Lysis Buffer plus 1 mM phenyl-methylsulfonyl fluoride (PMSF) to each plate (10 cm in diameter) and incubate the plate on ice for 5 minutes.
4. Scrape cells off the plate and transfer to an appropriate tube. Keep on ice.
5. Sonicate lysates on ice.
6. Microcentrifuge for 10 minutes at 4°C and transfer the supernatant to a new tube. The supernatant is the cell lysate. Store at –80°C in single-use aliquots.

C Test Procedure

1. After the microwell strips have reached room temperature, break off the required number of microwells. Place the microwells in the strip holder. Unused microwells must be resealed and stored at 4°C immediately.
2. Add 100 μl of Sample Diluent (supplied in each PathScan® Sandwich ELISA Kit, blue color) to a microcentrifuge tube. Transfer 100 μl of cell lysate into the tube and vortex for a few seconds. (Sample applied to the well can be diluted 1:1 when the suggested cell lysis buffer is used for cell extraction.) Individual data sheets for each kit provide information regarding an appropriate dilution factor for lysates and kit assay results.

3. Add 100 μl of each diluted cell lysate to the appropriate well. Seal with tape and press firmly onto top of microwells. Incubate the plate for 2 hours at 37°C. Alternatively, the plate can be incubated overnight at 4°C, which gives the best detection of target protein.
4. Gently remove the tape and wash wells:
 - a. Discard plate contents into a receptacle.
 - b. Wash 4 times with 1X Wash Buffer, 200 μl each time for each well.
 - c. For each wash, strike plates on fresh towels hard enough to remove the residual solution in each well, but do not allow wells to completely dry at any time. Clean the underside of all wells with a lint-free tissue.
5. Add 100 μl of Detection Antibody (green color) to each well. Seal with tape and incubate the plate for 1 hour at 37°C.
6. Repeat wash procedure as in Step 4C.
7. Add 100 μl of HRP-Linked secondary antibody (red color) to each well. Seal with tape and incubate the plate for 30 minutes at 37°C.
8. Repeat wash procedure as in Step 4C.
9. Add 100 μl of TMB Substrate to each well. Seal with tape and incubate the plate for 10 minutes at 37°C or 30 minutes at 25°C.
10. Add 100 μl of STOP Solution to each well. Shake gently for a few seconds.

NOTE: Initial color of positive reaction is blue, which changes to yellow upon addition of STOP Solution.

11. Read results.
 - a. Visual Determination — Read within 30 minutes after adding STOP Solution.
 - b. Spectrophotometric Determination — Wipe underside of wells with a lint-free tissue. Read absorbance at 450 nm within 30 minutes after adding STOP Solution.