

#2758 Store at -20°C

ISG15 (22D2) Rabbit mAb



✓ 100 µl
(10 western blots)

Orders ■ 877-616-CELL (2355)
orders@cellsignal.com
Support ■ 877-678-TECH (8324)
info@cellsignal.com
Web ■ www.cellsignal.com

rev. 09/09/10

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

Applications	Species Cross-Reactivity*	Molecular Wt.	Isotype
W, IP Endogenous	H, Hr	15 kDa	Rabbit IgG**


Background: Interferon-stimulated 15 kDa protein (ISG15), also known as ubiquitin cross-reactive protein (UCRP), is a member of the ubiquitin-like protein family and functions in various biological pathways from pregnancy to innate immune responses (1). Expression of ISG15 is stimulated by cellular exposure to type 1 interferons α and β , in addition to infection with viruses such as influenza B (2,3). After exposure to type I interferons, both lymphocytes and monocytes, in addition to some fibroblasts and epithelial cells, release ISG15 into culture medium (1,4). ISG15 has been shown to function as a cytokine, stimulating interferon γ secretion by monocytes and macrophages, proliferation of natural killer cells, and chemotactic responses in neutrophils (4,5). ISG15 has also been shown to function intracellularly, being covalently conjugated to other proteins by E1 (Ube1L), E2 (UbcH8) and E3 ligases via a multi-step process analogous to ubiquitination (6,7). ISG15 is removed from proteins by the ubiquitin processing protease Ubp43 (8). ISG15-protein conjugation (ISGylation) is induced by type 1 interferons, and target proteins include the serine protease inhibitor Serpin 2A, PLC γ 1, ERK1/2, Jak1 and Stat1 (9,10). Unlike ubiquitination, ISGylation does not target proteins for degradation, rather ISGylation increases Jak1 and Stat1 activity, enhancing the cellular response to interferons (11).

Specificity/Sensitivity: ISG15 (22D2) Rabbit mAb detects endogenous levels of the uncleaved precursor form of ISG15 protein. This antibody does not recognize the activated (cleaved) or conjugated forms of ISG15. The antibody does not cross-react with other ubiquitin family members, including ubiquitin, SUMO-1, SUMO-2, SUMO-3 and NEDD8.

Source/Purification: Polyclonal antibodies are produced by immunizing animals with a synthetic peptide corresponding to amino acids from human ISG15 protein. Antibodies are purified by protein A and peptide affinity chromatography.

Background References:

- Ritchie, K.J. and Zhang, D.E. (2004) *Semin. Cell Dev. Biol.* 15, 237–246.
- Korant, B.D. et al. (1984) *J. Biol. Chem.* 259, 14835–14839.
- Haas, A.L. et al. (1987) *J. Biol. Chem.* 262, 11315–11323.
- Knight, E. and Cordova, B. (1991) *J. Immunol.* 146, 2280–2284.
- D’Cunha, J. et al. (1996) *Proc. Natl. Acad. Sci. USA* 93, 211–215.
- Loeb, K.R. and Haas, A.L. (1992) *J. Biol. Chem.* 267, 7806–7813.
- Zhao, C. et al. (2005) *Proc. Natl. Acad. Sci. USA* 102, 10200–10205.
- Malakhov, M.P. et al. (2002) *J. Biol. Chem.* 277, 9976–9981.
- Malakhov, M.P. et al. (2003) *J. Biol. Chem.* 278, 16608–16613.
- Hamerman, J.A. et al. (2002) *J. Immunol.* 168, 2415–2423.
- Malakhova, O.A. et al. (2003) *Genes Dev.* 17, 455–460.



Western blot analysis of lysates from HeLa and A459 cells, untreated (-) or treated with IFN- α (1000 U/mL) (+) for 24 hours, using ISG15 (22D2) Rabbit mAb.

Entrez-Gene ID #9636
Swiss-Prot Acc. #P05161

Storage: Supplied in 10 mM sodium HEPES (pH 7.5), 150 mM NaCl, 100 µg/ml BSA, 50% glycerol and less than 0.02% sodium azide. Store at -20°C. Do not aliquot the antibody.

*Species cross-reactivity is determined by western blot.

**Anti-rabbit secondary antibodies must be used to detect this antibody.

Recommended Antibody Dilutions:

Western blotting	1:1000
Immunoprecipitation	1:25

For application specific protocols please see the web page for this product at www.cellsignal.com.

Please visit www.cellsignal.com for a complete listing of recommended companion products.

IMPORTANT: For western blots, incubate membrane with diluted antibody in 5% w/v BSA, 1X TBS, 0.1% Tween-20 at 4°C with gentle shaking, overnight.

Applications Key: W—Western IP—Immunoprecipitation IHC—Immunohistochemistry ChIP—Chromatin Immunoprecipitation IF—Immunofluorescence F—Flow cytometry E-P—ELISA-Peptide
Species Cross-Reactivity Key: H—human M—mouse R—rat Hm—hamster Mk—monkey Mi—mink C—chicken Dm—D. melanogaster X—Xenopus Z—zebrafish B—bovine
Dg—dog Pg—pig Sc—S. cerevisiae Ce—C. elegans Hr—horse All—all species expected Species enclosed in parentheses are predicted to react based on 100% homology.

© 2010 Cell Signaling Technology, Inc. Rabbit monoclonal antibody is produced under license (granting certain rights including those under U. S. Patents No. 5,675,063 and 7,429,487) from Epitomics, Inc.