

Rabbit (DA1E) mAb IgG XP™ Isotype Control (Alexa Fluor® 647 Conjugate)

100 µl

Orders ■ 877-616-CELL (2355)
orders@cellsignal.com

Support ■ 877-678-TECH (8324)
info@cellsignal.com

Web ■ www.cellsignal.com

rev. 01/27/11


This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

Concentration	Applications	Isotype
100 µg/ml	F	Rabbit IgG


Background: Isotype control antibodies are used to estimate the non-specific binding of target primary antibodies due to F_c receptor binding or other protein-protein interactions. An isotype control antibody should have the same immunoglobulin type as the test antibody.

Description: This Cell Signaling Technology antibody was conjugated to Alexa Fluor® 647 fluorescent dye and tested in-house for direct flow cytometric analysis of human cells.

Specificity/Sensitivity: Rabbit (DA1E) mAb IgG XP™ Isotype Control is not expected to recognize any specific antigen. It functions as an isotype control for rabbit monoclonal antibodies.



Flow cytometric analysis of Jurkat cells, using Rabbit (DA1E) mAb IgG XP™ Isotype Control (Alexa Fluor® 647 Conjugate) (red) compared to Phospho-S6 Ribosomal Protein (Ser235/236) (D57.2.2E) XP™ Rabbit mAb (Alexa Fluor® 647 Conjugate) #4851, untreated (blue) or IFN- α -treated (green).



Flow cytometric analysis of Jurkat cells, using XRabbit (DA1E) mAb IgG XP™ Isotype Control (Alexa Fluor® 647 Conjugate) (red) compared to Phospho-Akt (Ser473) (D9E) XP™ Rabbit mAb (Alexa Fluor® 647 Conjugate) #4075, untreated (green) or LY294002- and wortmannin-treated (blue).

Storage: Supplied in PBS (pH 7.2), less than 0.1% sodium azide, 2 mg/ml BSA. Store at 4°C. *Protect from light. Do not freeze.*

Note: This control antibody must be diluted to the same concentration (not dilution) as the specific antibody in analysis. See Directions on Use.

Recommended Antibody Dilutions:

Flow Cytometry:
Match IgG concentration of primary antibody.

Directions for Use: Important! This control antibody must be diluted to the same concentration (not dilution) as the specific antibody in analysis. Higher background fluorescence may result if excessive amounts of rabbit IgG isotype control are used. Do not use this antibody at 1:50 dilution in the same way as recommended for other CST Alexa Fluor® conjugated antibodies.

For application specific protocols please see the web page for this product at www.cellsignal.com.

Please visit www.cellsignal.com for a complete listing of recommended companion products.

The Alexa Fluor® dye antibody conjugates in this product are sold under license from Molecular Probes, Inc. for research use only, except for use in combination with DNA microarrays. The Alexa Fluor® dyes (except for Alexa Fluor® 430 dye) are covered by pending and issued patents.

Alexa Fluor® is a registered trademark of Molecular Probes, Inc.