

Ero1-L α Antibody

✓ 100 μ l
(10 western blots)

Orders ■ 877-616-CELL (2355)
orders@cellsignal.com
Support ■ 877-678-TECH (8324)
info@cellsignal.com
Web ■ www.cellsignal.com


rev. 03/16/10

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

Entrez-Gene ID #30001
Swiss-Prot Acc. #Q96HE7

Applications	Species Cross-Reactivity*	Molecular Wt.	Source
W Endogenous	H	60 kDa	Rabbit**

Background: Secretory proteins translocate into the endoplasmic reticulum (ER) after their synthesis where they are post-translationally modified and properly folded. To reach their native conformation, many secretory proteins require the formation of intra- or inter-molecular disulfide bonds (1). This process is called oxidative protein folding. Several oxidoreductases of the protein disulfide isomerase (PDI) family essential for disulfide formation and isomerization are localized to the ER (2). Studies have found that the ER-residing protein endoplasmic oxidoreductin-1 (Ero1) provides the oxidizing potential to the ER in *Saccharomyces cerevisiae* (3). *In vitro* experiments demonstrated that Ero1 is oxidized by molecular oxygen in a FAD-dependent manner and the oxidized Ero1 in turn serves as an oxidant for PDI (4). Two human homologs of Ero1, Ero1-like (Ero1-L α and β) have been identified (2,5). Ero1-L α is an ER membrane-associated N-glycoprotein that promotes oxidative protein folding and has been shown to be expressed in several cell lines and tissues (2).



Western blot analysis of extracts from various cell lines, using Ero1-L α Antibody.

Specificity/Sensitivity: Ero1-L α Antibody detects endogenous levels of total Ero1-L α protein.

Source/Purification: Polyclonal antibodies are produced by immunizing animals with a synthetic peptide corresponding to residues surrounding Leu218 of human Ero1-L α . Antibodies are purified by protein A and peptide affinity chromatography.

Background References:

- (1) Huppa, J.B. and Ploegh, H.L. (1998) *Cell* 92, 145–148.
- (2) Cabibbo, A. et al. (2000) *J. Biol. Chem.* 275, 4827–4833.
- (3) Frand, A.R. and Kaiser, C.A. (1998) *Mol. Cell* 1, 161–170.
- (4) Tu, B.P. et al. (2000) *Science* 290, 1571–1574.
- (5) Pagani, M. et al. (2000) *J. Biol. Chem.* 275, 23685–23692.

Storage: Supplied in 10mM HEPES (pH 7.5), 150 mM NaCl, 100 μ g/ml BSA and 50% glycerol. Store at -20°C. Do not aliquot the antibody.

*Species cross-reactivity is determined by western blot.


**Anti-rabbit secondary antibodies must be used to detect this antibody.

Recommended Antibody Dilutions:

Western blotting 1:1000

For application specific protocols please see the web page for this product at www.cellsignal.com.

Please visit www.cellsignal.com for a complete listing of recommended companion products.



IMPORTANT: For western blots, incubate membrane with diluted antibody in 5% w/v BSA, 1X TBS, 0.1% Tween-20 at 4°C with gentle shaking, overnight.

Applications Key: W—Western IP—Immunoprecipitation IHC—Immunohistochemistry ChIP—Chromatin Immunoprecipitation IF—Immunofluorescence F—Flow cytometry E-P—ELISA-Peptide
Species Cross-Reactivity Key: H—human M—mouse R—rat Hm—hamster Mk—monkey Mi—mink C—chicken Dm—D. melanogaster X—Xenopus Z—zebrafish B—bovine
 Dg—dog Pg—pig Sc—S. cerevisiae Ce—C. elegans Hr—Horse All—all species expected Species enclosed in parentheses are predicted to react based on 100% homology.