

Btk (C82B8) Rabbit mAb

✓ 100 µl
(10 western blots)

Orders ■ 877-616-CELL (2355)
orders@cellsignal.com
Support ■ 877-678-TECH (8324)
info@cellsignal.com
Web ■ www.cellsignal.com

rev. 11/11/11

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

Entrez-Gene ID #695
Swiss-Prot Acc. #Q06187

Applications	Species Cross-Reactivity*	Molecular Wt.	Isotype
W Endogenous	H, M	77 kDa	Rabbit IgG**


Background: Bruton's tyrosine kinase (Btk) is a member of the Btk/Tec family of cytoplasmic tyrosine kinases. Like other Btk family members, it contains a pleckstrin homology (PH) domain, and Src homology SH3 and SH2 domains. Btk plays an important role in B cell development (1,2). Activation of B cells by various ligands is accompanied by Btk membrane translocation mediated by its PH domain binding to phosphatidylinositol-3,4,5-trisphosphate (3-5). The membrane-located Btk is active and associated with transient phosphorylation of two tyrosine residues, Tyr551 and Tyr223. Tyr551 in the activation loop is transphosphorylated by the Src family tyrosine kinase, leading to autophosphorylation at Tyr223 within the SH3 domain, which is necessary for full activation (6,7). The activation of Btk is negatively regulated by PKCβ through phosphorylation of Btk at Ser180, which results in reduced membrane recruitment, transphosphorylation and subsequent activation (8). The PKC inhibitory signal is likely to be a key determinant of the B-cell receptor signaling threshold to maintain optimal Btk activity (8).

Specificity/Sensitivity: Btk (C82B8) Rabbit mAb detects endogenous levels of total Btk. This antibody does not cross-react with other unrelated proteins.

Source/Purification: Monoclonal antibody is produced by immunizing animals with a synthetic peptide corresponding to residues surrounding Pro80 of human Btk.

Background References:

- (1) Khan, W.N. (2001) *Immunol. Res.* 23, 147–156.
- (2) Lewis, C.M. et al. (2001) *Curr. Opin. Immunol.* 13, 317–325.
- (3) Salim, K. et al. (1996) *EMBO J.* 15, 6241–6250.
- (4) Rameh, L.E. et al. (1997) *J. Biol. Chem.* 272, 22059–22066.
- (5) Varnai, P. et al. (1999) *J. Biol. Chem.* 274, 10983–10989.
- (6) Rawlings, D.J. et al. (1996) *Science* 271, 822–825.
- (7) Park, H. et al. (1996) *Immunity* 4, 515–525.
- (8) Kang, S.W. et al. (2001) *EMBO J.* 20, 5692–5702.



Western blot analysis of extracts from various cell types using Btk (C82B8) Rabbit mAb.

Storage: Supplied in 10 mM sodium HEPES (pH 7.5), 150 mM NaCl, 100 µg/ml BSA and 50% glycerol. Store at -20°C. Do not aliquot the antibody.

*Species cross-reactivity is determined by western blot.

**Anti-rabbit secondary antibodies must be used to detect this antibody.

Recommended Antibody Dilutions:

Western blotting 1:1000

For application specific protocols please see the web page for this product at www.cellsignal.com.

Please visit www.cellsignal.com for a complete listing of recommended companion products.

IMPORTANT: For western blots, incubate membrane with diluted antibody in 5% w/v BSA, 1X TBS, 0.1% Tween-20 at 4°C with gentle shaking, overnight.

Rabbit monoclonal antibody is produced under license (granting certain rights including those under U. S. Patent No. 5,675,063) from Epitomics, Inc.