

Store at -20°C
#3999

Streptavidin-HRP

1 ml



Orders ■ 877-616-CELL (2355)
orders@cellsignal.com
Support ■ 877-678-TECH (8324)
info@cellsignal.com
Web ■ www.cellsignal.com

rev. 10/15/10

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

Applications	Species Cross-Reactivity
W	All

Description: This Cell Signaling Technology product is useful for the detection of biotinylated proteins (1,2). Conjugation of horseradish peroxidase (HRP) to streptavidin is obtained by cross linking the amino groups on streptavidin with the carbohydrate groups on HRP.

Background: Streptavidin is a 53 kDa homotetramer isolated from *Streptomyces avidinii* for use in isolation and detection bioassays (3). Each streptavidin subunit forms high affinity non-covalent bonds with the vitamin biotin. Because of its strong non-covalent interaction with biotin, streptavidin can be used to detect and isolate biotinylated proteins (1,2).


Chemiluminescent detection systems have emerged as the best all-around detection method for use with western blots and ELISA. These detection assays eliminate the hazards associated with radioactive materials and toxic chromogenic substrates. The speed and sensitivity of these methods are unequalled by traditional alternatives. Streptavidin-HRP is used with biotinylated proteins and specific chemiluminescent substrates to generate light signal. Streptavidin-HRP conjugates have a very high turnover rate, coupling high sensitivity with short reaction times

Specificity/Sensitivity: Streptavidin has a remarkably high affinity for its natural ligand, biotin. The complex and irregular structure of the biotin-binding site makes it highly optimized for biotin binding and confers great specificity to the streptavidin-biotin complexes (4).

Source/Purification: Streptavidin-HRP is purified by gel filtration. The conjugate product is free from unconjugated streptavidin and HRP.

Background References:

- (1) Updyke, T.V. and Nicolson, G.L. (1984) *J Immunol Methods* 73, 83-95.
- (2) Buckie, J.W. and Cook, G.M. (1986) *Anal Biochem* 156, 463-72.
- (3) Chaiet, L. and Wolf, F.J. (1964) *Arch Biochem Biophys* 106, 1-5.
- (4) Reznik, G.O. et al. (1998) *Proc Natl Acad Sci USA* 95, 13525-30.



Western blot analysis of extracts from HeLa cells, untreated or treated with Paclitaxel #9807, using Phospho-Histone H3 (Ser10) (D2C8) Rabbit mAb (Biotinylated) #3642 and Streptavidin-HRP #3999 for detection.

Swiss-Prot Acc. #P22629

Storage: Supplied in 136 mM NaCl, 2.6 mM KCl, 12 mM sodium phosphate (pH 7.4) dibasic, 2 mg/ml BSA, 50% glycerol. Store at -20°C. Do not aliquot antibody.

Recommended Antibody Dilutions:
Western blotting 1:2000

For application specific protocols please see the web page for this product at www.cellsignal.com.

Please visit www.cellsignal.com for a complete listing of recommended companion products.

IMPORTANT: Dilute Streptavidin-HRP in 5% nonfat dry milk, 1X TBS, 0.1% Tween 20. Incubate for 1 hour at room temperature shaking.

© 2010 Cell Signaling Technology, Inc.

Applications Key: W—Western IP—Immunoprecipitation IHC—Immunohistochemistry ChIP—Chromatin Immunoprecipitation IF—Immunofluorescence F—Flow cytometry E-P—ELISA-Peptide
Species Cross-Reactivity Key: H—human M—mouse R—rat Hm—hamster Mk—monkey Mi—mink C—chicken Dm—D. melanogaster X—Xenopus Z—zebrafish B—bovine
Dg—dog Pg—pig Sc—S. cerevisiae Ce—C. elegans Hr—Horse All—all species expected