

#4968 Store at -20°C

Pan-Actin Antibody



✓ 100 µl
(10 western blots)

Orders ■ 877-616-CELL (2355)
orders@cellsignal.com
Support ■ 877-678-TECH (8324)
info@cellsignal.com
Web ■ www.cellsignal.com

rev. 02/03/11

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.


Applications	Species Cross-Reactivity*	Molecular Wt.	Source
W, IHC-P Endogenous	H, M, R, Mk, Z (Pg, Dm, X, B)	45 kDa	Rabbit**

Background: Actin, a ubiquitous protein in eukaryotes, is the major component of the cytoskeleton. At least six isoforms are known in mammals. Nonmuscle β- and γ-actin, also known as cytoplasmic actin, are predominantly expressed in nonmuscle cells, controlling cell structure and motility (1). α-cardiac and α-skeletal actin are expressed in striated cardiac and skeletal muscles, respectively; two smooth muscle actins, α- and γ-actin, are found primarily in vascular smooth muscle and enteric smooth muscle, respectively. These actin isoforms regulate contractile potentials for muscle cells (1). Actin exists mainly as a fibrous polymer, F-actin. In response to cytoskeletal reorganizing signals during processes such as cytokinesis, endocytosis, or stress, cofilin promotes fragmentation and depolymerization of F-actin, resulting in an increase in the monomeric globular form, G-actin (2). The Arp2/3 complex stabilizes F-actin fragments and promotes formation of new actin filaments (2). It has been reported that actin is hyperphosphorylated in primary breast tumors (3). Cleavage of actin under apoptotic conditions has been observed *in vitro* and in cardiac and skeletal muscle (4-6). Actin cleavage by caspase-3 may accelerate ubiquitin/proteasome dependent muscle proteolysis (6).


Specificity/Sensitivity: Pan-Actin Antibody detects endogenous levels of total actin (all isoforms). The antibody also detects the 30 kDa actin fragment cleaved at glutamate 107.

Source/Purification: Polyclonal antibodies are produced by immunizing animals with a synthetic peptide corresponding to residues surrounding Asp244 of human β-actin. Antibodies are purified by protein A and peptide affinity chromatography.

- Background References:**
- (1) Herman, I.M. (1993) *Curr. Opin. Cell Biol.* 5, 48–55.
 - (2) Condeelis, J. (2001) *Trends Cell Biol.* 11, 288–293.
 - (3) Lim, Y.P. et al. (2004) *Clin. Cancer Res.* 10, 3980–3987.
 - (4) Kayalar, C. et al. (1996) *Proc. Natl. Acad. Sci. USA.* 93, 2234–2238.
 - (5) Communal, C. et al. (2002) *Proc. Natl. Acad. Sci. USA.* 99, 6252–6256.
 - (6) Du, J. et al. (2004) *J. Clin. Invest.* 113, 115–123.



Western blot analysis of extracts from HeLa, L929, C6 and COS cells using Pan-Actin Antibody.



Immunohistochemical analysis of paraffin-embedded human leiomyoma using Pan-Actin Antibody.

Entrez-Gene ID #60
Swiss-Prot Acc. #P60709

Storage: Supplied in 10 mM sodium HEPES (pH 7.5), 150 mM NaCl, 100 µg/ml BSA and 50% glycerol. Store at -20°C. Do not aliquot the antibody.

*Species cross-reactivity is determined by western blot.
**Anti-rabbit secondary antibodies must be used to detect this antibody.

Recommended Antibody Dilutions:
Western Blotting 1:1000
Immunohistochemistry (Paraffin) 1:50†
Unmasking buffer: Citrate
Antibody diluent: SignalStain® Antibody Diluent #8112
Detection reagent: SignalStain® Boost (HRP, Rabbit) #8114
†Optimal IHC dilutions determined using SignalStain® Boost IHC Detection Reagent.

For application specific protocols please see the web page for this product at www.cellsignal.com.
Please visit www.cellsignal.com for a complete listing of recommended companion products.

IMPORTANT: For western blots, incubate membrane with diluted antibody in 5% w/v BSA, 1X TBS, 0.1% Tween-20 at 4°C with gentle shaking, overnight.

© 2011 Cell Signaling Technology, Inc. Cell Signaling Technology® is a trademark of Cell Signaling Technology, Inc.

Applications Key: W—Western IP—Immunoprecipitation IHC—Immunohistochemistry ChIP—Chromatin Immunoprecipitation IF—Immunofluorescence F—Flow cytometry E-P—ELISA-Peptide
Species Cross-Reactivity Key: H—human M—mouse R—rat Hm—hamster Mk—monkey Mi—mink C—chicken Dm—D. melanogaster X—Xenopus Z—zebrafish B—bovine
Dg—dog Pg—pig Sc—S. cerevisiae Ce—C. elegans Hr—Horse All—all species expected Species enclosed in parentheses are predicted to react based on 100% homology.