

#7077

Store at -20°C

Anti-rat IgG, HRP-linked Antibody

✓ 1 ml

 Cell Signaling
TECHNOLOGY®

Orders ■ 877-616-CELL (2355)
orders@cellsignal.com

Support ■ 877-678-TECH (8324)
info@cellsignal.com

Web ■ www.cellsignal.com

rev. 11/09/10

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

Description: Affinity purified goat anti-rat IgG (H&L) antibody conjugated to horseradish peroxidase.


Applications: This product has been optimized for use as a secondary antibody in chemiluminescent western blotting applications.

Storage: Supplied in 10 mM sodium HEPES (pH 7.5), 150 mM NaCl, 2.0 mg/ml bovine serum albumin (BSA) and 50% glycerol. Store at -20°C.

Recommended Antibody Dilution:

1:1000–1:3000

Please visit www.cellsignal.com for a complete listing of recommended companion products.



After the primary antibody is bound to the target protein, a complex with HRP-linked secondary antibody is formed. The LumiGlo® is added and emits light during enzyme catalyzed decomposition.