

Anti-rat IgG, HRP-linked Antibody

1 ml

Orders ■ 877-616-CELL (2355)
orders@cellsignal.com

Support ■ 877-678-TECH (8324)
info@cellsignal.com

Web ■ www.cellsignal.com

New 04/06

This product is for *in vitro* research use only and is not intended for use in humans or animals.

Description: Affinity purified goat anti-rat IgG (H&L) antibody conjugated to horseradish peroxidase.

Applications: This product has been optimized for use as a secondary antibody in chemiluminescent Western blotting applications.

Storage: Supplied in 10 mM sodium HEPES (pH 7.5), 150 mM NaCl, 0.1 mg/ml bovine serum albumin (BSA) and 50% glycerol. Store at -20°C.

Recommended Antibody Dilution:


1:1000–1:3000

Companion Products:

Prestained Protein Marker, Broad Range (Premixed Format) #7720

Biotinylated Protein Ladder Detection Pack #7727

20X LumiGLO® Reagent and 20X Peroxide #7003



After the primary antibody is bound to the target protein, a complex with HRP-linked secondary antibody is formed. The LumiGLO® is added and emits light during enzyme catalyzed decomposition.

*Avoid repeated exposure to skin (see enclosed MSDS, consult our website or call to request the MSDS from CST).