



Confocal immunofluorescent images of C2C12 cells serum starved and treated with or without insulin as indicated and labeled with Phospho-Akt (Ser473) Antibody (red). Actin filaments have been labeled with fluorescein phalloidin.

Background References:

- (1) Franke, T.F. (1997) *Cell* 88, 435–437.
- (2) Burgering, B.T. and Coffey, P.J. (1995) *Nature* 376, 599–602.
- (3) Franke, T.F. et al. (1995) *Cell* 81, 727–736.
- (4) Alessi, D.R. et al. (1996) *EMBO J.* 15, 6541–6551.
- (5) Sarbassov, D.D. et al. (2005) *Science* 307, 1098–1101.
- (6) Jacinto, E. et al. (2006) *Cell* 127, 125–137.
- (7) Cardone, M.H. et al. (1998) *Science* 282, 1318–1321.
- (8) Brunet, A. et al. (1999) *Cell* 96, 857–868.
- (9) Zimmerman, S. et al. (1999) *Science* 286, 1741–1744.
- (10) Cantley, L.C. et al. (1999) *Proc. Natl. Acad. Sci. USA* 96, 4240–4245.
- (11) Vlahos, C. et al. (1994) *J. Biol. Chem.* 269, 5241–5248.
- (12) Hajdуч, E. et al. (2000) *FEBS Lett.* 492, 199–203.
- (13) Cross, D.A. et al. (1995) *Nature* 373, 785–789.
- (14) Diehl, J.A. et al. (1998) *Genes Dev.* 12, 3499–3511.
- (15) Gesbert, F. et al. (2000) *J. Biol. Chem.* 275, 39223–39230.
- (16) Zhou, B.P. et al. (2001) *Nat. Cell Biol.* 3, 245–252.
- (17) Nave, B.T. et al. (1999) *Biochem. J.* 344, 427–431.
- (18) Manning, B.D. et al. (2000) *Mol. Cell* 4, 648–657.
- (19) Manning, B.D. et al. (2002) *Mol. Cell* 10, 151–162.
- (20) Inoki, K. et al. (2002) *Nat. Cell Biol.* 4, 648–657.



Flow cytometric analysis of LNCaP cells, untreated (green) or LY294002-treated (blue), using Phospho-Akt (Ser473) Antibody compared to a nonspecific negative control antibody (red).



Western blot analysis of extracts from NIH/3T3 cells, untreated or treated with PDGF for the indicated times, using Phospho-Akt (Ser473) Antibody (upper) or Akt Antibody #9272 (lower).