

#9804 Store at -20°C

ATP

Concentration: 10 mM

1 ml



Cell Signaling
TECHNOLOGY®

Orders ■ 877-616-CELL (2355)
orders@cellsignal.com

Support ■ 877-678-TECH (8324)
info@cellsignal.com

Web ■ www.cellsignal.com

rev. 09/21/04

Description: Adenosine-5'-triphosphate (ATP) supplied as a 10 mM solution in doubly distilled water as a disodium salt.

Notes on Use: For *in vitro* kinase assays, we recommend a final ATP concentration of 200 μ M.

Companion Products:

Kinase Buffer (10X) #9802

Cell Lysis Buffer (10X) #9803