

HTScan[®] Tyrosine Kinase Buffer (4X)

 15 ml


Orders ■ 877-616-CELL (2355)
orders@cellsignal.com

Support ■ 877-678-TECH (8324)
info@cellsignal.com

Web ■ www.cellsignal.com

rev. 09/12/07

This product is for *in vitro* research use only and is not intended for use in humans or animals.

Description: HTScan[®] Tyrosine Kinase Buffer (TKB) can be used in tyrosine kinase activity assays.

Solutions and Reagents:

1X Kinase Buffer

60 mM HEPES

5 mM MgCl₂

5 mM MnCl₂

3 μM Na₃VO₄

Note: TKB must be supplemented with ATP (#9804) and 2.5 mM dithiothreitol (DTT).

Storage: Store at -20°C.

Companion Products:

Kinase Buffer (10X) #9802

Cell Lysis Buffer (10X) #9803

ATP (10 mM) #9804

Tyrosine Kinase Substrate Screening Kit #7450